

DecFill® LAB-PP Semi-Automatic Peristaltic Pump



DecFill® LAB-PP Peristaltic Machine can be used as stand alone or integrated into a fully automatic machine. All pump components are placed inside the machine, making it suitable to pump aggressive fluids without the risk of cross contamination. With a filling accuracy of up to +/- 0.5 % this machine provides a precise fill.

DecFill® LAB-PP Semi-Automatic Peristaltic Pump



Disposable Product Path Kit
Customized Sterile Tubing
with Connections

Technical Data	DecFill® LAB-PP	DecFill® LAB-PP XXL
Dosing Range	0.1 - 100ml	100ml - 20L
Electrical Supply	110/220V, 50/60Hz	



- ① HMI
- ② Pump Head

